

XRD machine

Turning on machine


Step 1: Switch on the AC main switch in the UPS room



Step 2: switch on the main power of the compressor from the wall and wait 4min for pressure generation



Step 3: switch on the main power in the MCB for the XRD machine and every switch



Step 4: press the power button on the machine to turn on the light press the button at the end of the row



Note: you can observe a light signifying the machine being turned on




Step 4: wait until KV and MA show 0 0 respectively, the XRD wings become horizontal



Step 5: Switch on the main switch on the MCB for the chiller by moving the nozzle clockwise.

Note: Make sure the water level is in between the range.





Step 7: Wait until the temperature reaches between 21 and 22 in the machine display (wait 1 min for stability of the value).



Step 8: Rotate the XRD machine key 90° clockwise

Note: Before turning on the PC and working on the software, make sure the XRD machine displays the KV and mA are 30 and 10 respectively.





Sample preparation

Step 7: Sample preparation and sample placement

Step one: check for all materials in the box i.e. sample holder, support board, spatula, brush, paper, and glass. And clean all using Acetone or an available cleaning agent, if you have to change the paper you can change it as well from the drawer at the back of the machine





paper finding area

Step two : Put the card board on the bottom metal support and close it with the metal closure



Step three : put your fine-powered sample in the middle portion of the sampe holder



Step four: Smooth your sample using glass road in a uniform way and remove excess powder beside the sample container on the side metal using a brush.



Catalysis
Science &
Technology



Sustainable
Energy & Fuels



RSC
Sustainability



te below to
e journals




Sponsor of
resentation Prize at
Engineering: Enabling Trans
toward Sustainable Future (CHEMTSF
IIT Roorkee
08-10 September 2022

Submit your best research to


<https://rsc.li/journals>

Step five: Hold it using two fingers like shown in the picture and clean the bottom part of the sample holder for trace powder remains.




Software part


Step 1: Turn on pc then open the software by the name Data Collector PC password is : **FIST**




Step 2: On the opening window header part click on Instrument and click on connect to connect with the machine




Step 3: In the stage select spinner for the showing window after clicking ok



Step 4: On the next window you have to do nothing so just click ok




Step 5: Close the auto opened small window



Step 6: Before the operation, we have to maintain KV and mA 45 and 40 respectively so wait till the machine displays so to do that we follow the following step.

Step one : Click on measure at header part and select program



Step two: In the coming window select RampUp-RT then click on open

Execute Program


Measurement type

All

Name	Measurement Type	Description	Creation Date	Created by	Modification Date
Ramp Down RT	General batch		28-Mar-24 3:16:42 PM	User	28-Mar-24 3:16:42
RampUp-FlatStage	General batch		17-Jun-22 3:46:39 PM	User	17-Jun-22 3:46:39
RampUp-RT	General batch		15-Jun-22 5:48:55 PM	User	15-Jun-22 5:48:55


Open Browse...

Close


 Use the drop-down list to filter on measurement type.
You can sort the available measurements by clicking on the column headers.

Step three : Click Ok on the coming window

Note : you will see the increase the KV and mA increasing to 45 and 40





Step four : After it reaches 40 and 45 you can close both rump-up windows (which are displayed as one and two in the picture below)



Step7: Before opening the door , click on the **unlock** doors button (and you will hear a vent release sound) and open the door by holding tight (be so gentle on the handling and don't try to open the door till you hear the vent sound), put your sample in the frame of the sample holder, and rotate your sample holder just to crosscheck if it is properly placed in the appropriate holding place of the machine and close it carefully till you hear the magnet sound for attachment.




Note: Before putting the sample make sure the stage is downward to carry your sample, if not double click on the sample stage option then as shown in the picture, **uncheck** the lift-up checkbox, and click apply then after the stage is down click ok.





Step 7: put your sample stage up for analysis

To do that click on the lift stage option which is shown in the picture




And you will get a window for lift up or down button , remember that if the lift up check box is checked it will move the lift up ward and if its unchecked and applied its for moving the lift down ward so for now we will use lift upward so check the lift up checkbox , then apply and then ok.



Step 8: Analysis (selecting the program)

Step one: Measure and select the program and on the next window click on Browse



Execute Program

Measurement type

All

Name	Measurement Type	Description	Creation Date	Created by	Modification Date
Ramp Down RT	General batch		28-Mar-24 3:16:42 PM	User	28-Mar-24 3:16:42 PM
RampUp-FlatStage	General batch		17-Jun-22 3:46:39 PM	User	17-Jun-22 3:46:39 PM
RampUp-RT	General batch		15-Jun-22 5:48:55 PM	User	15-Jun-22 5:48:55 PM

Open

Browse...

Close



Use the drop-down list to filter on measurement type.
You can sort the available measurements by clicking on the column headers.

Open

← → ↕ ⬆ This PC > OS (C:) > PANalytical > Data Collector > Programs >

Organize ▾ New folder

★ Quick access

Desktop

Downloads

Documents

Pictures

15.04.24

16.04.24

Raj

Screenshots

OneDrive

This PC

3D Objects

Desktop

Name

Date modified

Type

Size

3-axes Cradle

06-May-22 4:40 PM

File folder

5-axes Cradle

06-May-22 4:40 PM

File folder

Capillary Spinner

06-Jun-22 1:29 PM

File folder

Programmable XY with Programmable Z Platforms

06-Jun-22 1:29 PM

File folder

Raj

25-Oct-22 3:54 PM

File folder

Reflection-Transmission Spinner 3.0

06-Jun-22 1:29 PM

File folder

Sample Stage for

Date created: 15-Jun-22 5:06 PM

Size: 779 KB

SAXS-WAXS Stag

Folders: Standard Measurement

06-Jun-22 1:29 PM

File folder

Ramp Down RT.x

Files: GIXRD-BBHD-FlatStage-Pixel.xrdmp, ...

28-Mar-24 3:16 PM

XRDMP File

9 KB

RampUp-FlatStage.xrdmp

17-Jun-22 3:46 PM

XRDMP File


9 KB

RampUp-RT.xrdmp

15-Jun-22 5:48 PM


XRDMP File

9 KB



Step two : sample analysis by selecting the program and saving

Step 1 : click on measure and select program on from the header



Step 2: Then click on Browse and you will see another window opening

Execute Program

Measurement type

All

Name	Measurement Type	Description	Creation Date	Created by	Modification Date
Ramp Down RT	General batch		28-Mar-24 3:16:42 PM	User	28-Mar-24 3:16:42
RampUp-FlatStage	General batch		17-Jun-22 3:46:39 PM	User	17-Jun-22 3:46:39
RampUp-RT	General batch		15-Jun-22 5:48:55 PM	User	15-Jun-22 5:48:55

Open

Browse...

Close



Use the drop-down list to filter on measurement type.
You can sort the available measurements by clicking on the column headers.

Step 3: Then select 'Raj' and select the third program like you selected in program setup portion, then click on the folder icon and click ok.

Start


Program

Name: ...\\Powder-BBHD-RT-Spinning-Big Holder-Pixel.xrdmp

Type: Absolute scan

Description:

File

Name: Powder-BBHD-RT-Spinning-Big Holder-Pixel_1.xrdm 

Folder: C:\\XRD Data

Comment:

Sample

ID:

Name:

Prepared by:

Position

Diffractometer

2Theta (°):	5.0000	Phi (°):	0.0	X (mm):	
Omega Offset (°):	0.0000	Chi (°):		Y (mm):	
Omega (°):	2.5000	Phi Offset (°):	-19.2	Z (mm):	

Reflection

Unit cell: h k l:

OK Cancel Help

Step 4: Go to local disk D and select the folder name “**XRD DATA**” and create your own folder (if its first time) and select it or select if folder by your name exists and give you a sample name and save the name. after saving, It will direct you back to the previous window , then click ok to start the analysis.

Step 5: Once the analysis is completed close both the progress and graph showing windows and , double click on “ **sample stage reflection-transmission spinner 3.0** “ from the left button options for stage up/down action and , uncheck the “**lift up**” checkbox then click on apply, when the ok button is active , click on ok.

Step 6 : While trying to remove the sample, first press the button to unlock the doors , then once you hear the vent sound, you can hold the door opening handler gently and open the door to remove the sample.

Data conversion using Xpert hicrocore software

After completion of analysis , open expert hicroscore software from the pc and open the xrd result file on the xpert software, edit necessary editings like background subtraction and save it as per your interest eg. (asc,raw etc..)

Machine closing

Step one : remove sample holder

step two : make the sample holder in the machine to “lift up” position.

Step three : go to measure and select program then select “ramp down” then click on open then on the upcoming windows select ok

Step four : wait for KV and mA reading till it show 30 and 10 respectively.

Step four : disconnect the software by going to instrument tab and click on disconnect button and click ok on the upcoming window. And close the software

Step five : go to machine and rotate the key anticlockwise 90 degree to turn off X-ray you will see the x-ray light getting off from side lights of the machine

Step six : machine power off by off button by long pressing (3 to 4 seconds) .

Step seven : tun off the chiller by reversing back 90 degree anticlockwise

Step eight : turn of the compreser from main switch

Step nine : turn off computer

Step ten : turn off power switches from MCB (xrd, chiller , and computer)

File --- > open program --- > browse -> raj -> spinning --> big holder --> browse --> go to local disk D -> go to XRD data -> save accordingly -> once the program is saved the program will start automatically.